AGN ENTERPR















movement and high precision motorized specimen feed, for high throughput histology labs. Powerful sectioning motor enable us to cut even the hard specimen easily. The smartly located feature keypad enables us to control the operation easily. Apart from other safety feature it also has an emergency stop button. **Salient Features** • LCD Display to indicate section Alarm to indicate start/end of travel Motorized forward and backward at ● Adjustable specimen retraction
Finger protection guard Section Counter ● Motorized sectioning with adjustable speed • Spacious waste tray **Technical Specification** ● Section Thickness Setting Range: 0.5-100µm thickness, travel etc ● Trimming Thickness Setting Range:1-600μm ● Setting Values: 0.5-5μm increment step 0.5µm two different speeds ● 5µm-20μm increment step 1μm ● 20μm-60um increment step 5μm ● 60μm-100μm increment step 10μm ● Minimum setting value :0.5μm ● Sectioning Speed: Adjustable

Read More

SKU:

Price: ₹750.000.00

Stock: instock

Categories: Histopathology

Product Description

Incorporates a combination of motorized hand wheel movement and high precision motorized specimen feed, for high throughput histology labs. Powerful sectioning motor enable us to cut even the hard specimen easily. The smartly located feature keypad enables us to control the operation easily. Apart from other safety feature it also has an emergency stop button. Salient **Features** ● LCD Display to indicate section ● Alarm to indicate start/end of travel ● Motorized forward and backward at ● Adjustable specimen retraction ● Finger protection guard ● Section

AGN ENTERPRISES



MANUFACTURER & SUPPLIER
ALL KINDS OF LABORATORY EQUIPMENTS



100μm increment step 10μm → Minimum setting value :0.5μm → Sectioning Speed: Adjustable